

GUANGZHOU TOPMEDI CO., LTD

ENJOY LIFE WITH INTELLIGENCE

Seating position

✓ Standing position

Lying flat

Driving freely

Automatic Standing Wheelchair





Address: RM1505-1507, No. 83, Huadi Avenue Liwan Dist,. Guangzhou, 510380, China















Led Light

Safety Belt


Perfect Standing and Seating position





Lying







Technical Specification:

- Fully padded seat, back rest, neck rest and knee pads in easy to clean leather look fabric.
- Fully adjustable for various body shapes.
- Safely supports up to 120 kg/ 275.5 lbs
- Pneumatic tires (50 psi) for top comfort on all terrains
- 360 degree joystick controller
- Maximum distance on fully charged battery 50km/31 miles
- 3 x 3500N and 1 x 6000N motors for reliable power
- Slop capacity of up to 12 degree
- Fully balanced wheelchair, with intelligent electromagnetic brake (electronic regenerative disc brake)
- 4 x 12V / 20ah Lithium batteries-Chargeable or off board

Overall Length : 112cm/44inch
Overall Width : 66cm/26inch
Overall height : 110cm/43inch
Max load : 125kg/275.5lbs

Pneumatic tires : 50 psi

Front : 18 inch rim width

Back : 4 inch

Max slope capacity : 12 degree Seat depth : 45 cm/18 inch Seat width : 45 cm/ 18 inch Seat height : 56 cm/ 22 inch

Height of extended back rest : 60 cm/ 23.5 inch Arm rest height : 36 cm/ 14 inch

Net weight : 125 kg/ 276 lbs